
Rx ONLY Caution! Federal (USA) law restricts this device to sale by or on the order of a dentist.

CLEANING, DISINFECTING AND STERILIZING
NETTOYAGE, DÉSINFECTION ET STÉRILISATION
REINIGEN, DESINFIZIEREN UND STERILISIEREN
LIMPIEZA, DESINFECCIÓN Y ESTERILIZACIÓN
PULIZIA, DISINFEZIONE E STERILIZZAZIONE
RENGÖRING, DESINFICERING OCH STERILISERING

Reprocessing Instructions

DENTAL

EMS 
ELECTRO MEDICAL SYSTEMS

PLEASE READ THIS BEFORE BEGINNING WORK!

Please read these operating instructions carefully as they explain all the most important details and procedures. Please pay special attention to the safety precautions.

Always keep this instruction close at hand.

To prevent injury to people and damage to property, please heed the corresponding directives. They are marked as follows:

Danger
Risk of injury

Caution
Risk of property damage or environmental harm

Please note
Useful additional information and hints

Prohibited

Authorized

Must be sterilized at up to 135°C in the autoclave

Suitable for thermal disinfection

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Lämplig för värmedesinfektion

About this manual

The instructions in this manual are only applicable to the equipment it was delivered with.

Note that the English version of this manual is the master from which translations derive. In case of any discrepancy, the binding version is the English text.

We would be pleased to answer your questions or listen to your suggestions. We do, of course, provide support in case of technical problems. Please call directly your approved EMS repair center or your dealer.

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EMS

INTRODUCTION

These reprocessing instructions provide instructions for cleaning, disinfection, sterilization and packaging of EMS reusable products intended to be reprocessed in medical facilities.

The goal of reprocessing reusable products is to reduce bioburden and to achieve sterility of those products in order to eliminate the risk of product reuse related infection. Decisions regarding cleaning, disinfecting or sterilizing EMS medical and dental instruments are based on the potential risk of infection associated with their use.

It is recommended to use steam sterilization.

Remember that sterilization or high-level disinfection cannot be achieved unless the elements of the assembly are cleaned first.

If there is anything in the following instructions that is not clear, do not hesitate to contact EMS in Nyon.

If you find that the reprocessing instructions from the manufacturer seem to be inadequate, please inform EMS in Nyon about those inadequacies.

We encourage you to report adverse events related to device reprocessing. Report such events directly to EMS in Nyon.

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CONTENTS

Reprocessing instructions for reusable products	2
Application.....	2
Preparation.....	4
Basic principles	4
Post-operative treatment	6
Steps	6
Cleaning/disinfection: Manual and mechanical reprocessing.....	8
Mechanical cleaning and disinfection.....	10
Steps for mechanical cleaning and disinfection with a disinfectant	14
Manual procedure	16
Steps for manual cleaning	18
Steps for manual disinfection.....	20
Inspection and maintenance.....	22
Packaging and sterilization	24
Service life	30

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REPROCESSING INSTRUCTIONS FOR REUSABLE PRODUCTS

The instructions are binding for the reprocessing of all **reusable products*** of EMS. When necessary, additional product-specific instructions are included with the product to provide additional information.

*(Here after called "products")

Application

Important

Before use, carefully read the operating instructions of the EMS instrument and devices with which the product will be used.

Reusable products must be cleaned and sterilized prior to first use. They must be replaced after the number of operations specified by the manufacturer. Disposable products (single use) cannot be reused.

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PREPARATION

Basic principles

It is only possible to carry out effective sterilization after the completion of effective cleaning and disinfection. Please ensure that, as part of your responsibility for the sterility of products during use, only sufficiently validated equipment and product-specific procedures are used for cleaning/disinfection and sterilization, and that the validated parameters are adhered to during every cycle.

Please also observe the applicable legal requirements in your country as well as the hygiene regulations of the hospital or clinic. This applies especially with regard to the additional requirements for the inactivation of prions.

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POST-OPERATIVE TREATMENT

 The post-operative treatment must be carried out immediately, no later than 30 minutes after the completion of the operation. Additional information, where necessary, is provided in the respective product-specific usage instructions.

Steps

- 1 Rinse away any surface soiling on the product with distilled deionized water or with a cleaning agent.
- 2 Rinse through all lumina (e.g. irrigation and aspiration connection) at least 3 times in the normal direction of flow (no back rinsing) using a disposable syringe (min. volume 50 ml) filled with distilled/deionized water applied to the back nozzle.
- 3 An aldehyde-free cleaning and disinfection solution that is compatible with the products may also be used as an alternative rinsing solution. In such cases, it is necessary to rinse thoroughly afterwards at least 3 times with distilled, or deionized water.

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MECHANICAL CLEANING AND DISINFECTION

Thermal disinfection should be used if this function is available on your disinfectant. Use if possible a disinfecting cycle compliant with the standard EN ISO 15883.

Note that there is a risk of disinfectant residue on products when using chemical disinfectants.

Ensure the following criteria are met when selecting a disinfectant system:

- Disinfectant is proven effective through testing (e.g. FDA approved or CE marked / EN ISO 15883 compliant).
- Disinfectant has suitable baskets to hold small fragile products and has rinsing connections for the attachment to product lumina.
- The cleaning program is suitable for products to be processed and the rinsing cycle is sufficient.
- Only low microbe count (<10 cfu/ml) distilled or deionized water is used for all rinsing steps. (E.g. Aqua purificata, as per the specifications of Pharm. Eur. or USP).
- Air used for drying is HEPA filtered.
- Disinfectant is serviced and checked on a regular basis.

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Ensure the following criteria are met when selecting a cleaning and disinfection agent:

- Chemicals are compatible with products.
- With non-thermal disinfection, a suitable disinfectant with tested effectiveness (e.g. FDA approved or CE marked), that is compatible with the cleaning agent, must be used.

 Concentrations and contact times specified by the manufacturer of the cleaning and disinfection agent must be followed. Only freshly prepared solutions may be used.

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Steps for mechanical cleaning and disinfection with a disinfectant

- 1 Carefully place the products in the disinfection basket. Fastening of the products is only permissible if they are freely moveable in the fixture. The products are not permitted to make contact with one another.
- 2 Using a suitable rinsing adaptor, connect the product lumina to the rinsing connections of the disinfectant.
- 3 Start the program.
- 4 Remove the products from the disinfectant and start the inspection (see section Inspection and maintenance) after the program ends.
- 5 Wrap the products directly following disinfection and drying (see section Packaging and sterilization). If necessary, repeat drying of the product in a clean place.

Verification of the fundamental suitability of the products for effective cleaning and disinfection by machine was provided by an independent, accredited testing laboratory when used with the G 7736 CD disinfectant, Miele & Cie. GmbH & Co., Gutersloh, (thermal disinfection), and the cleaning agent Neodisher® Medizym (Dr. Weigert GmbH & Co. KG, Hamburg). The procedure described above was adhered to as part of this.

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MANUAL PROCEDURE

When selecting the cleaning and/or disinfection agent to be used, ensure that:

- These are fundamentally suitable for the cleaning and/or disinfection of the products and compatible with one another,
- The disinfection agent holds a tested effectiveness (e.g. DGHM, FDA approval or CE marking),
- The chemicals used are compatible with the products.

It is absolutely essential that the concentrations and contact times specified by the manufacturer of the cleaning and disinfection agent are adhered to. Only freshly prepared solutions may be used. The disinfectant solution is not permitted to foam.

Only sterilized or low microbe count distilled/deionized water (< 10 cfu/ml) can be used for all rinsing steps. Also ensure the presence of a sufficiently low endotoxin and particle concentration (e.g. Aqua purificata, as per the specifications of Pharm. Eur. or USP).

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Steps for manual cleaning

- 1 Completely disassemble the handpiece and instruments, if applicable.
- 2 Place the products in the cleaning solution for at least the minimum the time and the concentration specified by the manufacturer of the cleaning/disinfection agent.
- 3 Remove any externally-attached soiling by brushing carefully with a soft brush or a soft cloth.
- 4 Rinse the products vigorously at least five times, each time with fresh distilled or deionized water (each product lumen with at least 50 ml of water). Repeat the cleaning process if the last rinsing does not run clear, or if stains are still visible on the product.

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Steps for manual disinfection

- 1 Place the products in the disinfection solution for at least the minimum the time specified by the manufacturer of the disinfection agent. Use an aldehyde-free disinfection solution for cleaning general and surgical instruments.
- 2 If applicable: rinse all lumens of the instrument five times at the beginning and at the end of the soaking time with a single use syringe (minimal volume of 50 ml).
- 3 Remove the instrument from the cleaning solution and post rinse them intensively with a low contaminated and deionized water
- 4 Dry the products with filtered compressed air (max. 3 bars).
- 5 Wrap the products directly following disinfection and drying (see section Packaging and sterilization). If necessary, repeat drying of the product in a clean place.

Verification of the fundamental suitability of the products for effective manual cleaning and disinfection was provided by an independent, accredited testing laboratory for use of the following cleaning/disinfection agents:

For handpiece FD 350 (from Dürr GmbH) and for the instruments ID 212 (from Dürr GmbH). The procedure described above was adhered to as part of this.

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INSPECTION AND MAINTENANCE

 If stains are still visible on the product after cleaning/disinfection, the entire cleaning/disinfection procedure must be repeated. Products with visible damage, chip/flake loss, corrosion or bent out of shape must be disposed of (no further use is permissible).

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PACKAGING AND STERILIZATION

Do not exceed the maximum number of sterilization cycles.

Only cleaned and disinfected products are permitted to be sterilized. Prior to sterilization, the products need to be placed in a suitable sterilization container:

- Compliant with EN ISO 11607,
- Resistant to 138°C, with adequate steam permeability,
- Maintained on a regular basis.

If double, single-use sterilization packaging (double bag) is to be used, this must also comply with EN ISO 11607 and be suitable for steam sterilization (temperature resistant to 138°C with adequate steam permeability).

Use only the following listed steam sterilization procedures for sterilization; other sterilization procedures are not permissible:

- Fractional pre-vacuum procedure*

* *Fractional pre-vacuum procedure = steam sterilization with repetitive pre-vacuum.*

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- Steam sterilizer in accordance with EN 13060 or EN 285 validated in compliance with EN ISO 17665,
- Maximum sterilization temperature 138°C,
- Sterilization time at least 3 min. at 132/134°C (fractional pre-vacuum procedure),
- Sterilization at 134°C for a maximum of 20 minutes is permissible.

Verification of the fundamental suitability of the products for effective steam sterilization was provided by an independent, accredited testing laboratory for use of the fractionated vacuum procedure.

(EuroSelectomat from MMM Münchener Medizin Mechanik GmbH, Planegg, validated according to DIN EN 285). The procedure described above was adhered to as part of this.

The hot-air sterilization and radio-sterilization procedure may not be used (as it causes the destruction of products).

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The manufacturer assumes no responsibility for the use of other sterilization procedures (e.g. ethylene oxide, formaldehyde and low temperature plasma sterilization). In such cases, please observe the respective valid standards (EN ISO 14937/ANSI AAMI ISO 14937 or the procedure-specific standard) and verify the suitability and effectiveness in principle of the procedure (if necessary, including investigations on sterilizing agent residue), taking into account the specific product geometry as part of the validation.

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SERVICE LIFE

The products have been designed for a large number of sterilization cycles. The materials used in their manufacture were selected accordingly. However with every renewed preparation for use, thermal and chemical stresses will result in ageing of the products.

If the number of permissible re-sterilization cycles is restricted, this will be pointed out in the product specific instructions.

The use of ultrasound baths and strong cleaning and disinfection fluids (alkalizing pH > 9 or acid pH < 5) can reduce the life span of products. The manufacturer accepts no liability in such cases.

The products may not be exposed to temperatures above 138°C.

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EMS Corporation